

## Packaging options for all your shipping requirements.

For your convenience Labcabs can supply IATA approved packaging for Biological substances CAT B, and Infectious substances CAT A, plus suitable packaging for shipping frozen, chilled or ambient goods.

The following examples can be purchased from Labcabs as part of the total freight booking process.

### Bio Bottle



Bio Bottle sizes: 500ml or 2 litre  
Packaging instruction 650 for shipping UN3373 samples.  
Packaging instruction 602 for shipping UN2814 samples.

### Medium Cool Cube Shipper Box



10 Litre Polyurethane container designed to ship products that are temperature sensitive.  
Ideal for maintaining temperatures from ( chilled) +2 C to +8 C or (frozen) -78 C for long transit periods, up to 96 hours.

### Cool 5 Shipper Box



5 Litre Polyurethane container designed to ship products that are temperature sensitive.  
Ideal for maintaining temperatures from ( chilled) +2 C to +8 C or (frozen) -78 C for long transit periods, up to 48 hours

### Ambient Shipper Box & Bio Pouch




Suitable for Domestic and International shipping at room temperature when combined.

Bio Pouch pressure tested to 95kPa

Medium: 24cm x 34cm

Large: 43cm x 34cm

### Chilled Food Container



Quality Polystyrene container ideal for shipping food samples, chilled or frozen. Domestic or International.

Dry ice top ups in transit included in Labcabs service.

### Data Temperature Logger



Log Tag is a simple to use electronic device that records temperatures until the item is delivered.  
The unit has a unique external temperature sensor for fast reaction time. Rapid download for reporting details, easy to use, reliable, credit card size, alert indicator, inspection mark, high performance.

**Call 1300 522 222 for all Bookings and queries on your next shipment.**

Labcabs Pty Ltd PO Box 276 Cannon Hill QLD 4170

[www.labcabs.com.au](http://www.labcabs.com.au)